

TSE

Smart technologies



BABY BLUE LIGHT FTL – 600

Phototherapy lamp

Baby Blue Light FTL – 600 Phototherapy lamp

The phototherapy lamp Baby Blue Light FTL – 600 represents the most common method of adequate treatment of neonatal hyperbilirubinemia. It brings an economical solution of adequate treatment along with a perfect price/performance ratio.

The simple suspension of the lamp on a height-adjustable stand enables both tilting up to 90° in one direction, and effective irradiation of the patient's entire body on an infant bed or a warming bed. Anti slip pads on the lamp surface ensure safe location on an incubator cover after removing the lamp from the stand.

A counter with a digital display records the number of operation hours and thus ensures safe and efficient operation.

Power voltage	230 V ± 10 %, 50 Hz
Max. Input power of the unit	170 VA
Input power of light sources	2 x 20 W
Light source	discharge tube PHILIPS (blue)
Wave length	400 – 550 nm
Light source life	approx. 2 000 hours

Dimensions	
Height	115 mm
Length	820 mm
Width	210 mm

Weight	4.5 kg
---------------	--------

Optional accessories	
Stand FTL – 600	adjustable height 1025 – 1475 mm



TSE spol. s r. o.
Mánesova 74, 371 52 České Budějovice, Czech Republic
T: + 420 386 721 149, F: + 420 386 721 102
medical@tse.cz, www.tse.cz

